

1. Part Number: A406-3000

2. Photograph:



3. Product Information:

Description	Straight Sided Jug 3000ml
Material of Construction:	316L stainless steel
Nominal Volume:	3000ml
Brimful Volume:	3350ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	223mm
Top Diameter	147mm
Top Diameter – including spout:	161mm
Max. Width – from spout to handle	194mm
Nominal Weight:	1080g
Graduation Marks?	Yes (every 500ml starting from 1000ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Oct-20