

**1. Part Number:** A406-5000

**2. Photograph:**



### 3. Product Information:

Description	Straight Sided Jug 5000ml
Material of Construction:	316L stainless steel
Nominal Volume:	5000ml
Brimful Volume:	5470ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	245mm
Top Diameter	178mm
Top Diameter – including spout:	192mm
Max. Width – from spout to handle	225mm
Nominal Weight:	1760g
Graduation Marks?	Yes (every 500ml starting from 1000ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 19-Apr-22