


**1. Part Number:** A416-10

**2. Photograph:**



Dimensions

A = 119mm

B = 98mm

C = 363mm

D = 212mm

Wall Thickness  
= 1.2mm

### 3. Product Information:

Description:	Stainless Container 10L
Materials of Construction:	- Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade)
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	10 Litres
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Weight:	3.6 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Aug-20



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.  
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com