


**1. Part Number:** A416-10

**2. Photograph:**



Dimensions

A = 119mm

B = 98mm

C = 363mm

D = 212mm

Wall Thickness  
= 1.2mm

### 3. Product Information:

Description:

Stainless Container 10L

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)

Method of Construction:

All welds ground and polished  
Crevice free interior

Surface Finish:

Better than 0.5 microns Ra

Nominal Volume

10L

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Weight:

3.6 Kg (body + lid + gasket + clamp)

Recommended Storage Conditions:

Dry and ambient temperature

Pressure Rating:

None – these containers are not pressure rated vessels

Notes:

- The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3
- These containers are primarily designed for storing powders
- They are not suitable for the storage or transportation of beer

## 4. Spare Gaskets:

Supplied with a 4" silicone gasket as standard. Gaskets are consumable parts and will need to be inspected and changed at the first signs or wear.

Part No.	Description
K201-4	Gasket 4" silicone – FDA acceptable grade
K202-4	Gasket 4" EPDM – FDA acceptable grade
K203-4	Gasket 4" Viton – FDA acceptable grade

## 5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Apr-26