


1. Part Number: A416-15

2. Photograph:



3. Product Information:

Description:	Stainless Container 15L
Materials of Construction:	<ul style="list-style-type: none"> - Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade)
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	15L
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Weight:	4.6 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Pressure Rating:	None – these containers are not pressure rated vessels
Notes:	<ul style="list-style-type: none"> - The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3 - These containers are primarily designed for storing powders - They are not suitable for the storage or transportation of beer

4. Spare Gaskets:

Supplied with a 4" silicone gasket as standard. Gaskets are consumable parts and will need to be inspected and changed at the first signs or wear.

Part No.	Description
K201-4	Gasket 4" silicone – FDA acceptable grade
K202-4	Gasket 4" EPDM – FDA acceptable grade
K203-4	Gasket 4" Viton – FDA acceptable grade

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Apr-26