


1. Part Number: A416-25

2. Photograph:



Dimensions

A = 119mm

B = 98mm

C = 460mm

D = 302mm

Wall Thickness
= 1.2mm

3. Product Information:

Description:	Stainless Container 25L
Materials of Construction:	- Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade)
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	25 Litres
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Weight:	5.3 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 31-Mar-21



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com