


**1. Part Number:** A416-30

**2. Photograph:**



Dimensions

A = 119mm

B = 98mm

C = 602mm

D = 302mm

Wall Thickness  
= 1.2mm

### 3. Product Information:

Description:

Stainless Container 30L

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)

Method of Construction:

All welds ground and polished  
Crevice free interior

Surface Finish:

Better than 0.5 microns Ra

Nominal Volume

30 Litres

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Weight:

6.9 Kg (body + lid + gasket + clamp)

Recommended Storage Conditions:

Dry and ambient temperature

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Aug-20



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.  
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com