


1. Part Number: A417-5

2. Photograph:



3. Product Information:

Description	Stainless Container 5L with External GL45 Thread
Materials of Construction:	<ul style="list-style-type: none"> - Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade) - Cap: 316L stainless steel - Cap Seal: Silicone (FDA acceptable grade)
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	5 Litres
Brimful Volume:	5.7 Litres
Dimensions	A = 119mm B = 98mm C = 316mm D = 164mm
Weight:	2.7 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Note:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Oct-20