


1. Part Number: A418-15

2. Photograph:



3. Product Information:

Description: Stainless Container 15L with Internal GL45 Thread

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)
- Cap: 316L stainless steel
- Cap Seal: Silicone (FDA acceptable grade)

Method of Construction: All welds ground and polished
Crevice free interior

Surface Finish: Better than 0.5 microns Ra

Nominal Volume: 15 L

Brimful Volume: 17.5 L

Dimensions

- A = 119mm
- B = 98mm
- C = 360mm
- D = 275mm

Weight: 4.60 Kg (body + lid + gasket + clamp)

Recommended Storage Conditions: Dry and ambient temperature

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25