


**1. Part Number:** A418-1

**2. Photograph:**



### 3. Product Information:

Description	Stainless Container 1L with Internal GL45 Thread
Materials of Construction:	<ul style="list-style-type: none"> <li>- Body: 316L stainless steel</li> <li>- Lid: 316L stainless steel</li> <li>- Clamp: 304 stainless steel</li> <li>- Gasket: Silicone (FDA acceptable grade)</li> <li>- Cap: 316L stainless steel</li> <li>- Cap Seal: Silicone (FDA acceptable grade)</li> </ul>
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	1 Litre
Dimensions	<ul style="list-style-type: none"> <li>A = 119mm</li> <li>B = 101mm</li> <li>C = 152mm</li> <li>D = 104mm</li> </ul>
Weight:	1.6 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Note:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25