


1. Part Number: A418-2

2. Photograph:



3. Product Information:

Description	Stainless Container 2L with Internal GL45 Thread
Materials of Construction:	<ul style="list-style-type: none"> - Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade) - Cap: 316L stainless steel - Cap Seal: Silicone (FDA acceptable grade)
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	2 Litres
Dimensions	<ul style="list-style-type: none"> A = 119mm B = 98mm C = 227mm D = 124mm
Weight:	2.0 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Note:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20