


**1. Part Number:** A418-5

**2. Photograph:**



**3. Product Information:**

Description	Stainless Container 5L with Internal GL45 Thread
Materials of Construction:	<ul style="list-style-type: none"> <li>- Body: 316L stainless steel</li> <li>- Lid: 316L stainless steel</li> <li>- Clamp: 304 stainless steel</li> <li>- Gasket: Silicone (FDA acceptable grade)</li> <li>- Cap: 316L stainless steel</li> <li>- Cap Seal: Silicone (FDA acceptable grade)</li> </ul>
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	5 Litres
Brimful Volume:	5.7 Litres
Dimensions	A = 119mm B = 98mm C = 316mm D = 164mm
Weight:	2.7 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Note:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25