


1. Part Number: A418-8

2. Photograph:



3. Product Information:

Description	Stainless Container 8L with Internal GL45 Thread
Materials of Construction:	<ul style="list-style-type: none"> - Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade) - Cap: 316L stainless steel - Cap Seal: Silicone (FDA acceptable grade)
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	8 litres
Brimful Volume:	8.1 litres
Dimensions	A = 119 mm B = 98 mm C = 283 mm D = 212 mm
Weight:	3.15 Kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Note:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
 Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com