

**1. Part Number:** A426-1000**2. Photograph:****3. Product Information:**

Description:	Mini-Pot 1000ml
Method of Construction:	Crevice free
Material of Construction:	316L Stainless Steel
Nominal Volume:	1000ml
Brimful Volume:	1000ml
Overall Height (body & lid):	154mm
Body Height:	150mm
Maximum Body Diameter:	100mm
Lid Diameter:	102mm
Nominal Weight:	1535g
Finish:	Bright Polish (Ra < 0.8 microns)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: If you require container with a watertight seal please see our containers type A836

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25