

**1. Part Number:** A426-250

**2. Photograph:**



**3. Product Information:**

Description:	Mini-Pot 250ml
Method of Construction:	Crevice free
Material of Construction:	316L Stainless Steel
Nominal Volume:	250ml
Brimful Volume:	250ml
Overall Height (body & lid):	73mm
Body Height:	71mm
Maximum Body Diameter:	74mm
Lid Diameter:	77mm
Nominal Weight:	530g
Finish	Bright Polish (Ra < 0.8 microns)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: If you require container with a watertight seal please see our containers type A836

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25