

**1. Part Number:** A426-7

**2. Photograph:**



### 3. Product Information:

Description:	Mini-Pot 7ml
Method of Construction:	Crevice free machined from solid bar (no welds)
Material of Construction:	316L Stainless Steel
Nominal Volume:	7ml
Overall Height (body & lid):	45mm
Body Height:	42mm
Maximum Body Diameter:	16mm
Lid Diameter:	19mm
Nominal Weight:	23g
Finish	Machined Finish (Ra < 1.0 microns)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25