

1. Part Number: A436-750

2. Photo:



3. Product Information:

Description	Filling Scoop 750ml
Nominal Volume:	750ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	313mm
Body Length:	183mm
Body Diameter:	101mm
Body Metal Thickness:	1.5mm
Handle Length:	130mm
Handle Outer Diameter:	30mm
Handle Interior Diameter:	26mm
Nominal Weight:	660g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Jan-25