


**1. Part Number:** A446-15

**2. Photograph:**



A = 285 mm  
B = 275 mm  
C = 305 mm  
D = 267 mm  
E = 275 mm

### 3. Product Information:

Description:	Mini-Drum 15L
Method of Construction:	Body welds ground and polished Patch plates for swing handle and lid handle are spot welded
Material of Construction:	316L Stainless Steel (body and lid)
Nominal Volume:	15L
Nominal Weight:	3.3 kg (body & lid)
Surface Finish:	Bright Polish (Ra < 0.5 microns)
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 20-Nov-20