


1. Part Number: A456-200-26

2. Photo:



3. Product Information:

Description	Funnel 200mm top diameter, 26mm spout diameter
Use	Suitable for low viscosity liquids
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Dimensions	A - Top Diameter 200mm B - Diameter of Spout 26mm C - Length of Spout 35mm D - Overall Height 190mm E - Wall Thickness 1.5mm to 2mm
Internal Angle:	60°
Nominal Weight:	720g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Nov-25