


	<h1>Product Specification</h1>	Specification Number	A466-50
		Revision	01

Part No. A466-50

Description Storage Container 50L

Materials of Construction: - Body: 316L stainless steel
- Lid: 316L stainless steel



Nominal Volume: 50 litres

Method of Construction: Body welds ground and polished
Side Handles are spot welded to body
Patch plate for lid handle is spot welded to lid


Surface Finish: Better than 0.5 microns Ra

Dimensions
A = 560 mm
B = 535 mm
C = 435 mm
D = 350 mm
E = 355 mm

Weight: 7.45 Kg (body + lid)

Graduation Marks: Yes (every 10L)

Handle Type: Side Handles

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 2-Jul-14
----------------------------	--	---