


	<h1>Product Specification</h1>	Specification Number	A466-75
		Revision	01

Part No. A466-75

Description Storage Container 75L

Materials of Construction: - Body: 316L stainless steel  
- Lid: 316L stainless steel



Nominal Volume: 75 litres

Method of Construction: Side Handles are spot welded to body  
Patch plate for lid handle is spot welded to lid


Surface Finish: Better than 0.5 microns Ra

Dimensions  
A = 625 mm  
B = 600 mm  
C = 480 mm  
D = 400 mm  
E = 405 mm

Weight: 9.8 Kg (body + lid)

Graduation Marks: Yes (every 10L)

Handle Type: Side Handles

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 17-Nov-16
Page 1 of 1	<b>Sampling Systems Ltd.</b> 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com	