

**1. Part Number:** A504-4

**2. Photograph:**



**3. Product Information:**

Description	Stainless Spoon 10ml 305mm
Material of Construction:	304 stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	10 ml
Overall Length:	305 mm
Bowl Length:	65 mm
Max. Bowl Width:	32 mm
Weight:	105 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Jan-25