

**1. Part Number:** A504-5

**2. Photograph:**



### 3. Product Information:

Description	Stainless Spoon 10ml 500mm
Material of Construction:	304 stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	10 ml
Overall Length:	500 mm
Bowl Length:	65 mm
Max. Bowl Width:	32 mm
Weight:	180 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Jan-25