

**1. Part Number:** A506-1000

**2. Photograph:**



### 3. Product Information:

Description	PharmaSpoon 10ml 1000mm
Material of Construction:	316L stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	10 ml
Overall Length:	1000 mm
Bowl Length:	70 mm
Max. Bowl Width:	35 mm
Weight:	371 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

### 4. Approval:

Signed:



Name: Kuldeep Flora  
 Position: Quality Manager  
 Date: 5-Dec-24