

1. Part Number: A506-210

2. Photograph:



3. Product Information:

| | |
|---------------------------|---|
| Description | PharmaSpoon 2ml 210mm |
| Material of Construction: | 316L stainless steel |
| Method of Construction: | Fully welded, the weld is ground and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Nominal Volume: | 2 ml |
| Overall Length: | 210 mm |
| Bowl Length: | 35 mm |
| Max. Bowl Width: | 20 mm |
| Weight: | 26 g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this scoop are of vegetable origin |
| Storage Conditions: | Dry, ambient conditions |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 5-Dec-24