

1. Part Number: A506-300

2. Photograph:



3. Product Information:

Description	PharmaSpoon 10ml 300mm
Material of Construction:	316L stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	10 ml
Overall Length:	300 mm
Bowl Length:	70 mm
Max. Bowl Width:	35 mm
Weight:	104 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 5-Dec-24