

1. Part Number: A566-02

2. Photograph:



3. Product Information:

Description:	TC Container 2L
Materials of Construction:	<ul style="list-style-type: none"> - Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade)
Method of Construction:	All welds ground and polished Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume	2000ml
Overall Height:	195mm
Body Height:	185mm
Body Diameter:	130mm
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Weight:	2.1kg (body + lid + gasket + clamp)
Recommended Storage Conditions:	Dry and ambient temperature
Notes:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3 These containers are primarily designed for storing powders. They are not suitable for the storage or transportation of beer.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Sep-24