

1. Part Numbers: A586-10

2. Photograph:



3. Product Information:

Description	Heavy Duty Bucket 10L
Nominal Volume:	10L
Brimful Volume:	10L
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free interior Lugs for handle are spot welded to exterior of bucket
Surface Finish:	Brush Finish - better than 1.0 microns Ra
Max. Diameter:	298mm
Bucket Height:	233mm
Graduations:	No
Skirt Band:	No
Nominal Weight:	2.1 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20