

**1. Part Numbers:** A586-15

**2. Photograph:**



**3. Product Information:**

Description	Heavy Duty Bucket 15L
Nominal Volume:	15L
Brimful Volume:	15L
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free interior Lugs for handle are spot welded to exterior of bucket
Surface Finish:	Brush Finish - better than 1.0 microns Ra
Max. Diameter:	330mm
Bucket Height:	305mm
Graduations:	No
Skirt Band:	No
Nominal Weight:	2.8 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20