

**1. Part Number:** A596-400

**2. Photograph:**



**3. Product Information:**

Description	Straight Tweezers 316Lss - 400mm
Material of Construction:	316L stainless steel
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	405mm
Blade Length	375mm
Nominal Weight:	205g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Mar-23