

1. Part Number: A606-5000

2. Photograph:



3. Product Information:

| | |
|---------------------------------|---|
| Description | GMP PharmaScoop® 5000ml |
| Nominal Volume: | 5000ml |
| Material of Construction: | 316L stainless steel |
| Method of Construction: | Crevice free, all welds ground and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Overall Length: | 500mm |
| Body Length: | 350mm |
| Body Diameter: | 195mm |
| Nominal Weight: | 1425g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this scoop are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25