

1. Part Number: A611-150

2. Photograph:



3. Product Information:

Description:	Laboratory Spoon 150mm
Material of Construction:	304 Stainless steel
Overall Length:	150mm
Bowl Length:	35mm
Bowl Width:	28mm
Nominal Bowl Volume:	2ml
Thickness of metal:	1.2mm
Nominal Weight:	21g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Jan-25