


1. Part Number: A611-210

2. Photograph:



3. Product Information:

Description:	Laboratory Spoon 210mm
Material of Construction:	304 Stainless steel
Overall Length:	212mm
Bowl Length:	40mm
Bowl Width:	30mm
Nominal Bowl Volume:	3ml
Thickness of metal:	1.2mm
Nominal Weight:	30g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Jan-25