

1. Part Number: A611-250

2. Photograph:



3. Product Information:

Description:	Laboratory Spoon 250mm
Material of Construction:	304 Stainless steel
Overall Length:	253mm
Bowl Length:	50mm
Bowl Width:	37mm
Nominal Bowl Volume:	6ml
Thickness of metal:	1.5mm
Nominal Weight:	56g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25