

1. Part Number: A611-300

2. Photograph:



3. Product Information:

Description:	Laboratory Spoon 300mm
Material of Construction:	304 Stainless steel
Overall Length:	303mm
Bowl Length:	55mm
Bowl Width:	43mm
Nominal Bowl Volume:	12ml
Thickness of metal:	1.5mm
Nominal Weight:	81g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Jan-25