

1. Part Number: A613-250

2. Photograph:



3. Product Information:

Description:	Flat Laboratory Spoon 250mm
Material of Construction:	304 Stainless steel
Overall Length:	255mm
Max. Width of Ends:	32mm
Thickness of metal:	1.5mm
Nominal Weight:	55g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Jan-25