

1. Part Number: A613-300

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description: | Flat Laboratory Spoon 300mm |
| Material of Construction: | 304 Stainless steel |
| Overall Length: | 306mm |
| Max. Width of Ends: | 40mm |
| Thickness of metal: | 1.5mm |
| Nominal Weight: | 79g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin. |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25