

1. Part Number: A618-235

2. Photograph:



3. Product Information:

Description:	Narrow Spoon Spatula 235mm
Material of Construction:	304 Stainless steel
Overall Length:	235mm
Bowl Length:	40mm
Bowl Width:	18mm
Nominal Bowl Volume:	2ml
Width of Flat End:	9mm
Length of Flat End:	45mm
Nominal Weight:	22g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Nov-25