


1. Part Number: A621-260

2. Photograph:



3. Product Information:

Description:	Beaker Tongs 260mm
Material of Construction:	304 Stainless steel
Overall Length:	260mm
Diameter of Arms:	6mm
Nominal Weight:	165g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Jul-20