


**1. Part Number:** A622-310

**2. Photograph:**



**3. Product Information:**

Description:	Non-Slip Beaker Tongs 310mm
Material of Construction:	- 304 Stainless steel - PVC coating on end of arms
Overall Length:	310mm
Diameter of Arms:	6mm
Max. Height of Arms:	22mm
Nominal Weight:	219g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for use with beakers over 45mm in diameter

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Jul-20