


1. Part Number: A626-8

2. Photograph:



A = 285 mm
B = 270 mm
C = 235 mm
D = 200 mm
E = 205 mm

3. Product Information:

Description:	Mini-Drum with Toggles - 8L
Method of Construction:	- Body welds ground and polished - Patch plates for toggles, swing handle and lid handle are spot welded
Materials of Construction:	- Body: 316L stainless steel - Lid: 316L stainless steel - Toggle Clamp: 304 stainless steel - O-ring: Silicone Rubber (conforms to FDA 177.2600)
Nominal Volume:	8L
Surface Finish:	Bright Polish (Ra < 0.5 microns)
Graduation Marks:	None
Handle Type:	Swing Handle
Number of Toggles	3
Nominal Weight:	2.5 kg (body + lid + toggles + O-ring)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25