

**1. Part Number:** A627-210

**2. Photograph:**



**3. Product Information:**

Description:	Chattaway Spatula 210mm
Material of Construction:	304 Stainless steel
Overall Length:	212mm
Max Width:	9mm
Blade Thickness:	1.0mm
Nominal Weight:	16g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 7-Jan-25