

1. Part Number: A627-250

2. Photograph:



3. Product Information:

Description: Chattaway Spatula 250mm
Material of Construction: 304 Stainless steel
Overall Length: 250mm
Max Width: 10mm
Blade Thickness: 1.0mm
Nominal Weight: 20g
BSE/TSE statement: All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25