

1. Part Number: A627-300

2. Photograph:



3. Product Information:

Description:	Chattaway Spatula 300mm
Material of Construction:	304 Stainless steel
Overall Length:	303mm
Max Width:	10mm
Blade Thickness:	1.0mm
Nominal Weight:	21g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25