

**1. Part Number:** A628-100

**2. Photograph:**



**3. Product Information:**

Description:	Micro Chattaway Spatula 100mm
Material of Construction:	304 Stainless steel
Overall Length:	103mm
Max Width:	4mm
Blade Thickness:	0.5mm
Nominal Weight:	2g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Jan-25