

1. Part Number: A630-120

2. Photograph:



3. Product Information:

Description:	Curved Spatula 120mm
Material of Construction:	304 Stainless steel
Overall Length:	120mm
Blade Length:	31mm
Blade Width:	15mm
Material Thickness:	1.2mm
Nominal Weight:	14g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25