


1. Part Number: A630-210

2. Photograph:



3. Product Information:

Description:	Curved Spatula 210mm
Material of Construction:	304 Stainless steel
Overall Length:	210mm
Blade Length:	51mm
Blade Width:	25mm
Material Thickness:	1.2mm
Nominal Weight:	39g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 8-Jan-25