

**1. Part Number:** A638-140

**2. Photograph:**



**3. Product Information:**

Description:	Straight Needle 140mm
Material of Construction:	304 Stainless steel
Overall Length:	147mm
Handle Width:	6mm
Needle Length:	40mm
Material Thickness:	2mm
Nominal Weight:	9g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product has a sharp needle point and should be used with caution

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 8-Jan-25