

1. Part Number: A640-140

2. Photograph:



3. Product Information:

Description:	Curved Needle 140mm
Material of Construction:	304 Stainless steel
Overall Length:	145mm
Handle Width:	6mm
Needle Length:	40mm
Material Thickness:	2mm
Nominal Weight:	9g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product has a sharp needle point and should be used with caution

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 8-Jan-25