

1. Part Number: A643-200

2. Photograph:



3. Product Information:

Description:	Industrial Scoop
Nominal Scoop Volume:	200ml
Material of Construction:	Stainless Steel
Method of Construction:	Fully welded handle
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	225mm
Body Length:	140mm
Max. Body Width:	85mm
Max. Body Height:	55mm
Nominal Weight:	155g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 17-May-22

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.